



Illustrations contain optional extras.

PARAVAN

PR HEAVY DUTY

Built to last

The strongest member in our range of power wheelchair solutions is the PR heavy duty. This power wheelchair model has been specially developed for maximum capacity. Weight limits up to 270 kg are no problem for the PR heavy duty. The chassis has been made wider, the wheel base lengthened, and the construction strengthened. Seats measuring up to 70 cm in width can now be fitted on the chassis. According to the latest studies, patients are getting heavier and heavier. We are catering for this development with the PR heavy duty.

YOUR BENEFITS



Highlight: Seat tilt module with stand-up aid:

Electric continuous seat tilt with stand-up aid, tilt angle up to 50° backwards and 20° forwards

Highlight: Electric footrest unit module:

Electric lower-leg angle adjustment up to 90°
Electric lower-leg length adjustment up to 120 mm
Individual footrests, spring mounted and foldable

Series model PR heavy duty:

- Front drive electric wheelchair
- Highly maneuverable thanks to width of only 72 cm
- Ground clearance approx. 70 mm
- Anti-tipping device located on the front
- Manual seat-length adjustment
- Individual footrests, spring mounted and foldable
- Puncture-resistant tyres with a grooved treadfront 3.00 x 8 inches and rear 2.00 x 50 dual-tyres with swing axle
- 10 A outdoor charger without fan
- Swing away control panel holder
- High performance and high-torque 550 W motors
- 2 x maintenance-free 88 Ah 12 V batteries
- Standard speed 4 mph
- Transport tie down points 2 x front and 1 x rear
- Maximum user weight: 270 kg (weight of person)

Key facts and options for PR heavy duty

- Made in Germany: the highest level of manufacturing and product quality. Developed by PARAVAN entirely in-house
- Tested for safety and maximum stability thanks to an extended wheelbase and wheel spacer
- Reinforced chassis
- Use of high-strength steel (electrogalvanised) used, preventing twisting
- Seat height manually adjustable
- Tilt function can be adjusted manually or electrically (option)
- Fixed, double-reinforced armrests for greater stability
- Perfect stability when driving thanks to dual tires
- Functional, convenient joystick control
- Numerous options for special constructions and accessories
- Diverse electrical functions available as options
- Customised modern paint finishes
- Vehicle docking station to secure the wheelchair in a vehicle (optional)