

PARAVAN PR 30/II

Our all-round talent

Discover the PARAVAN PR 30/II! The advanced wheelchair is modular, equipped with extensive function options and can be modified. The high-performance motor technology supports a maximum user weight of 140 kg. What's more, the PR 30 has a standard lifting function with a lifting height of up to 30 cm. The side supports of the innovative seat design have flexible width and angle adjustment. The side guards of the seat are movable, optionally removable and available in various heights. The bold chassis shape in combination with the lighting concept let the PR 30/II shine in a modern design.



IHRE VORTEILE



Highlight: Additional seat tilt module with stand-up aid:

Electric continuous seat tilt with stand-up aid, tilt angle up to 50° backwards

Highlight: Elektrische Fußstützeneinheit:

Electric lower-leg angle adjustment up to 90° Electric lower-leg length adjustment up to 120 mm Individual footrests, spring mounted and foldable

Basismodell PR 30/II:

- Vorderradgetriebener Elektrorollstuhl
- Nutzergewicht max. 140 kg (optional 200 kg)
- Zentrale Fußstützeneinheit mit einzelnen Fußplatten
- PARAVAN Kontur-Sitzsystem K620 mit neun individuell wählbaren Sitzgrößen (Lehne und Sitzteil), mit gesticktem PARAVAN-Schriftzug und Sitzkissen
- · Manual seat-length adjustment
- The armrests can be folded away and their height and inclination is adjustable; the head rest is also adjustable.
- Integrated back plate with thread inserts (e.g. for drip stand, pushing aid)
- Electric seat height adjustment (lifting function) lifting height 360 mm to 800 mm, infinite
- Modern lighting design and indicators advisable for side
- Joystick control with up to 6 functions and a swing arm mounting
- Puncture-proof tyres with grooved tread
- Very maneuverable, only 640 mm wide
- Ground clearance approx. 70 mm
- 10 A outdoor charger without fan
- Standard speed 4 mph
- Integrated spring mounted anti-tip protection
- Black tyres (non-marking)
- Rear wheels with suspension

Keyfacts und Optionen PR 30/II:

- Made in Germany: Top manufacturing and product quality. Developed entirely by PARAVAN
- Seat and backrest expansion (optional)
- Side bolsters on the seat part with adjustable width and angle (optional)
- Removable side bolsters in various heights (optional)
- Electrical backrest adjustment (optional)
- Electrical seat length adjustment (optional)
- Seat tilt 45° (optional)
- Modular system, can be adapted individually to meet the needs of the relevant clinical picture
- Large choice of individual seat pads (optional)
- Individual electric footrests with length adjustment (optional)
- Electric track gyro stabilisation (optional)
- Other controls available as an option (eg. chin control)
- Special controls & environment controls optional
- Distinctive modern colours
- Vehicle docking station to secure the wheelchair in a vehicle (optional)